

TECHNICAL SHEET

1. Product Identification

| | |
|---------------------|-------------------------------|
| Designation: | Toise Murale a Rouleau Rapide |
| SKU: | MU-TOISE-11 |
| Brand: | MediUnit Clinical |
| Category: | Diagnostics |

2. Technical Specifications

| | |
|-------------------------|------------|
| Packaging: | Unit |
| Plage De Mesure: | 0 - 220 cm |

3. Compliance & Documentation

| | |
|------------------------|--|
| CE Certificate: | Available on request or in the compliance pack |
| Sterilization: | Non-sterile / N.A. |
| Packaging: | Unit |

4. Clinical Description

Toise de mesure murale simple et efficace. Ruban metallique deroulant avec curseur de lecture precis.