

1. Product Identification

Designation:	Pese-Personne Medical Mecanique a Cadran
SKU:	MU-SCALE-09
Brand:	MediUnit Clinical
Category:	Diagnostics

2. Technical Specifications

Packaging:	Unit
Capacite Max:	150 kg
PreCision:	500g

3. Compliance & Documentation

CE Certificate:	Available on request or in the compliance pack
Sterilization:	Non-sterile / N.A.
Packaging:	Unit

4. Clinical Description

Pese-personne medical mecanique a cadran pour controle ponderal de routine en consultation, depistage ou suivi nutritionnel. - Le grand cadran et la base antiderapante facilitent une lecture rapide sur poste de passage.