

1. Product Identification

Designation:	Seringue Luer Slip 3 Pieces (Multi-Tailles)
SKU:	NEO176550SLIP
Brand:	Neomedic
Category:	Syringes

2. Technical Specifications

Packaging:	Unit
Volume:	50ml
Needle Gauge:	21G
Needle Length:	1.5
Connection Type:	Luer Slip
Material:	Polypropylene/Polyethylene
Sterility:	true
Single Use:	true
Latex Free:	true

3. Compliance & Documentation

CE Certificate:	Available on request or in the compliance pack
Sterilization:	Sterile EO
Packaging:	Unit

4. Clinical Description

La seringue Luer Slip 3 Pieces (Multi-Tailles), reference NEO176550SLIP, est conçue pour l'administration précise de fluides et de solutions cliniques. Sa conception en trois pièces améliore la maîtrise du dosage, la transparence du cylindre facilite le contrôle visuel et la connexion Luer Slip assure une liaison rapide avec les accessoires compatibles. Fabriquée en polypropylène de qualité médicale, stérilisée à l'oxyde d'éthylène (EO) et conditionnée individuellement, elle est destinée à un usage unique professionnel en clinique et à l'hôpital.